



Kavin Scientific Products

Price List 2016-2017



Glasswares



Liquid Handling



Protective Gears



Kavin Scientific Products

Certificate of Registration

This is to Certify that
Quality Management System of


KAVIN SCIENTIFIC PRODUCTS

7A, NATARAJAN STREET (SAIDAI SAMBANDAM ST),
SAIDAPET, CHENNAI - 600 015, TAMILNADU, INDIA

has been assessed and found to conform to the requirements of
ISO 9001:2008
for the following scope :

MANUFACTURE AND SUPPLIER OF LABORATORY
GLASSWARE AND SCIENTIFIC INSTRUMENTS

Certification No	1015QCE78		
Initial Registration Date	: 10/08/2015	Issuance Date	: 10/08/2015
Date of Expiry*	: 09/08/2018		
1st Surve. Due	: 10/07/2016	2nd Surve. Due	: 10/07/2017



Director



M4570910IN



JAS-ANZ



www.jas-anz.org/register

Absolute Quality Certification Pvt. Ltd.

Accredited By: Joint Accreditation System of Australia & New Zealand (JAS-ANZ)

401, Ashok Bhavan, 93, Nehru Place, New Delhi - 110 019, | +91. 11. 4105 1139 | e-mail : info@absolutequalitycertification.com.

*Validity of the Certificate is subject to successful completion of surveillance audit on or before of due date. (In case surveillance audit is not allowed to be conducted, this certificate shall be suspended/withdrawn).

Certificate Verification: Please Re-check the validity of certificate at <http://www.absolutequalitycertification.com/activeclients.aspx> or www.absolutequalitycertification.com at Active Clients.

Certificate is the property of Absolute Quality Certification Pvt. Ltd. and shall be returned immediately when demanded

ISO 9001 : 2008



Product Range:

Autoclaves	BOD Incubators	Blood Storage Cabinets	BioSafety Cabinet
Carbon & Sulphur Determination Apparatus	Constant Temperature bath	Cryostat Bath	Bunsen burners
Environmental Chamber	Frozen/Semen Storage cabinets	Fume Exhaust Hoods	Deep Freezers (Chest & Upright)
Furnace (Grooved)	Furnace Tubular (1200°c)	Furnace High Temperature (1400°c)	Flocculator (Jar Test Apparatus)
Humidity and Temp. Cabinet	Hot Plates	Heating Mantles	Hot Air Oven
Ice Flaking Machine	Incinerator	Incubator	Humidifier
Incubator Shaker	Inoculation hood	Kjeldahl (Digestion & Distillation)	Industrial Tray Drier
Low Temperature Bath	Mortuary Chamber	Muffle Furnace	Laminar Air Flow
Melting Point Apparatus	Orbital Shakers	Oil Baths	Magnetic Stirrers
Plant Growth Chamber	Rortary Shakers	Sieve Shakers	Oven(Industrial & Universal)
Sand baths	Tray Driers	Tissue Floatation bath	Soxhlet Extraction Units
VDRL Rotator	Vacuum Pump (Single./Double/Oil free)	Waterbaths	Vacuum Oven
Waterbath incubator shakers	Wrist action shakers	Water Still (wall & Table)	Waterbaths(Serological)



Manufactured & Marketed by:

Kavin Scientific Products



7A, Natarajan St., (Saidai Sambandan St.,)
 Saidapet, Chennai-600015, INDIA
 Tel: 044-43302975, 24714329
 E-Mail: kavinscientific@yahoo.co.uk,
 kavinscientific@gmail.com
 Website: www.kavinscientific.in



K-LAB Interio



Product categories:

- Wall bench
- Instrument table
- Anti Vibration Table
- Island bench
- Wall Storage cupboard
- Reagent Rack
- Chemical Storage Cupboard
- Bench type Fume Hood
- Distillation Fume Hood
- Walk in type fume hood.

CHENNAI

Kavin Scientific Products

An ISO9001-2008 Certified Company

7A, Natarajan St, (Saidai Sambandam St.)
Saidapet, Chennai-15.

Phone : 044-24714329 Telefax : 044-43302975

Cell : 98400 54224, 99400 94224

Email : kavinscientific@yahoo.co.uk

kavinscientific@gmail.com

MADURAI

233, South Marat Street,

South Gate, Madurai - 625001.

Cell : 98404 54224, 88838 38889

Email : kavinscientific.mdu@gmail.com

Website : www.kavinscientific.in



Kavin Scientific Products

PRICE LIST

HEAD OFFICE	BRANCH OFFICE
<p>CHENNAI 7A, Natarajan St, (Saidai Sambadhan St,) Saidapet, Chennai-15. Phone : 24714329, Tele Fax : 044-43302975 Mobile : 98400 54224, 99400 94224 Email : kavinscientific@yahoo.co.uk</p>	<p>MADURAI 233, South Marat Street, South Gate, Madurai-625001 Mobile : 98404 54224, 88838 38889 Email : kavinscientific.mdu@gmail.com</p>



LABORATORY GLASSWARES

(MADE FROM BOROSILICATE GLASS)



This type of borosilicate glass contains approx. 81% silicon dioxide and approx. 13% boron trioxide is featured in high chemical durability, low coefficient of linear thermal expansion and therefore even to high resistance against temperature changes

Borosilicate glass represent an internationally defined glass type, suitable for manufacturing the majority of laboratory glass, devices and industrial apparatus, The borosilicate glass features : Excellent corrosion resistance, smooth unporous surface, perfect transparency, catalytically indifferent, odourless and tasteless,

CHEMICAL DURABILITY

Borosilicate glass resists practically all aggressive agent except Hydrofluoric acid, Phosphoric acid and Concentrated Hot alkalis which considerably affect the surfaces of the glass. Hydrofluoric acid damages glass surface mostly at low concentrations and it is thus applied to etching scales on volumetric glass. Phosphoric acid and alkalis damage glass only slightly at low temperatures and concentrations. Continual changing aciduous and alkalinous medium increases corrosions.

REVIEW OF THE MOST IMPORTANT PHYSICAL CHARACTERISTICS

Medium linear coefficient of thermal expansivity α	20-300° C	3,3 10 ⁻⁶ K ⁻¹
Water resistance	98° C	class 1
Softening point (Littleton's)	10 ^{7.6} PAS	825°C
Density	20°C	2.23 g cm ⁻³


INDEX

ADAPTERS		DISHES	
Cone	- 2	Petri	- 16
Claisen Heads Sloping	- 6	Crystallizing	- 16
Expansion	- 1	Evaporating	- 16
Intermediate Receiver	- 5	EXTRACTION APPARATUS	
Multiple Two Necks	- 1	Soxhlet, Complete	- 17
Plain Still heads	- 5	Spare Allihn	- 17
Recovery bend	- 5	Spare Extractor	- 17
Reduction	- 1	FLASKS	
Receiver	- 3-4	Conical Narrow Mouth	- 17
Swan Neck	- 2	Conical Wide Mouth	- 17
Splash Heads	- 6	Conical with i/c Stopper	- 18
Steam Distillation Heads	- 6	Conical with Screwcap	- 18
With T Connection	- 2	Culture, Haffkine	- 21
ARSENIC DETERMINATION	- 33	Distilling	- 21
ASSEMBLIES DISTILLATION	- 29-34	Filter	- 21
BEAKERS		Flat Bottom, Plain	- 18
Low Form	- 8	Flat Bottom, Centre Neck	- 19
Tall Form	- 8	Iodine with i/c Stopper	- 18
BOTTLES		Kjeldahl	- 21-22
BOD	- 8	Pear Shape	- 22
Reagent	- 8-9	Pear Shape, for Rotary	- 23
Screwcap	- 9	Round Bottom, Plain	- 18
Specific Gravity	- 10	Round Bottom, Centre Neck	- 19-20
Wash	- 10	Round Bottom, Two Necks	- 20
Weighing tall form	- 10	Round Bottom, Three Necks	- 20
Weighing Squat form	- 10	Volumetric	- 22
BURETTES		FUNNELS	
Automatic with Reservoir	- 11	Buchner	- 25
Glass Stopcock	- 10	Long Stem	- 23
Glass Stopcock, Schellbach	- 11	Pressure Equalising	- 24
Rotoflow Teflon Screwcock	- 11	Short Stem	- 23
Rotoflow Teflon Schellbach	- 11	Separating, Pear	- 24
Teflon Stopcock	- 11	LIQUID HANDLING PRODUCTS	
CONDENSERS		Micropipettes, Fixed	- 35
Air	- 11	Micropipettes, Variable	- 35
Allihn	- 12	Pipette Controller	- 35
Double Surface	- 12	PIPETTES	- 25-26
Friedrichs	- 12	PROTECTIVE GEARS	
Graham Coil	- 12	Eye Protection	- 36
Columns, Chromotography	- 14-16	Eye Wash Units	- 38
Crow Receivers	- 14	Fire Protection	- 37
Crucibles, Filter	- 25	Hand Protection	- 37
CYLINDERS		Rspirator Protection	- 38
Nessler	- 13	PYKNO METER	- 8
Measuring	- 13	STOPPERS	- 28
DESSICATORS		Thermometer, Pocket	- 7
Plain	- 16	TILT MEASURE	- 26
Vacuum	- 16	TUBES	
		Air/Steam	- 7
		Adapters for Vacuum Release	- 7
		Adapters for Vacuum Control	- 7
		Blood Sugar	- 27
		Centrifuge	- 27
		Culture	- 27
		Drying	- 7
		Filter	- 26
		Test Tubes	- 28



0100 - Reduction Adapters,
Made from Heat Resistant Borosilicate Glass, KAVIN

Socket	Cone	Qty/pack	Price/Pc
14/23	19/26	5	140.00
14/23	24/29	5	168.00
14/23	29/32	5	200.00
19/26	24/29	5	184.00
19/26	29/32	5	218.00
19/26	34/35	5	235.00
19/26	45/40	5	295.00
19/26	55/44	5	335.00
24/29	29/32	5	235.00
24/29	34/35	5	250.00
24/29	45/40	5	315.00
24/29	55/44	5	330.00
29/32	34/35	5	265.00
29/32	45/40	5	290.00
29/32	55/44	5	500.00
34/35	45/40	5	335.00
34/35	55/44	5	500.00
45/40	55/44	5	545.00




0110 - Expansion Adapters,
Made from Heat Resistant Borosilicate Glass, KAVIN

Socket	Cone	Qty/pack	Price/Pc
19/26	14/23	5	140.00
24/29	19/26	5	200.00
24/29	14/23	5	168.00
29/32	24/29	5	235.00
29/32	19/26	5	218.00
29/32	14/23	5	218.00
34/35	29/32	5	265.00
34/35	24/29	5	300.00
24/35	19/26	5	300.00
34/35	14/23	5	300.00
45/40	34/35	5	450.00
45/40	29/32	5	450.00
45/40	24/29	5	450.00
55/44	45/40	5	630.00
55/44	34/35	5	630.00




0120 - Multiple Adapters, Two Necks Parallel,
Made from Heat Resistant Borosilicate Glass, KAVIN

Socket	Cone	Qty/pack	Price/Pc
14/23	14/23	4	210.00
19/26	19/26	4	268.00
19/26	24/29	4	275.00
19/26	34/35	4	350.00
24/29	24/29	4	335.00
29/32	29/32	4	430.00



0130 - Multiple Adapters, Two Necks, One Vertical and One At Angle
Made from Heat Resistant Borosilicate Glass, KAVIN

Socket	Cone	Qty/pack	Price/Pc
19/26	19/26	2	290.00
19/26	24/29	2	278.00





0140 - Multiple Adapters, Three Necks, Two Parallel and One At Angle

Made from Heat Resistant Borosilicate Glass, KAVIN



Parallel Neck Socket	Angle Neck Socket	Cone	Qty/pack	Price/Pc
14/23	14/23	19/26	2	335.00
19/26	14/23	24/29	2	385.00
19/26	19/26	24/29	2	385.00
19/26	19/26	29/32	2	395.00
19/26	19/26	34/35	2	470.00
19/26	19/26	55/44	2	660.00

0150 - Swan Neck Adapters, For use with Airleak Tube,


Made from Heat Resistant Borosilicate Glass, KAVIN



Socket	Cone	Qty/pack	Price/Pc
14/23	14/23	4	185.00
19/26	19/26	4	220.00
19/26	24/29	4	270.00
24/29	24/29	4	300.00

0160 - Adapters with 'T' Connection,

Made from Heat Resistant Borosilicate Glass, KAVIN



Socket	Cone	Qty/pack	Price/Pc
14/23	14/23	2	168.00
14/23	19/26	2	200.00
14/23	24/29	2	215.00
19/26	19/26	2	227.00
19/26	24/29	2	227.00
19/26	34/35	2	268.00

0170 - Cone Adapters Connection Plain,


Made from Heat Resistant Borosilicate Glass, KAVIN



Socket	Qty/pack	Price/Pc
10/19	4	76.00
14/23	4	84.00
19/26	4	100.00
24/29	4	109.00
29/32	4	126.00
34/35	4	175.00

0180 - Cone Adapters, Straight Connection with Stopcock,

Made from Heat Resistant Borosilicate Glass, KAVIN



Cone	Qty/pack	Price/Pc
14/23	4	250.00
19/26	4	260.00
24/29	4	294.00
29/32	4	300.00

0190 - Cone Adapters, Right Angle Connection Plain,

Made from Heat Resistant Borosilicate Glass, KAVIN



Cone	Qty/pack	Price/Pc
10/19	4	76.00
14/23	4	84.00
19/26	4	100.00
24/29	4	109.00
29/32	4	126.00
34/35	4	175.00




0200 - Cone Adapters, Right Angle Connection Stopcock,
 Made from Heat Resistant Borosilicate Glass, KAVIN

Cone	Qty/pack	Price/Pc
19/26	4	260.00
24/29	4	290.00
29/32	4	300.00
34/35	4	395.00



0210 - Socket Adapters, Straight Connection with Stopcock,
 Made from Heat Resistant Borosilicate Glass, KAVIN

Cone	Qty/pack	Price/Pc
14/23	4	200.00
19/26	4	260.00
24/29	4	290.00
29/32	4	395.00




0220 - Receiver Adapters, Plain Bend, Short,
 Made from Heat Resistant Borosilicate Glass, KAVIN

Socket	Qty/pack	Price/Pc
10/19	5	93.00
14/23	5	100.00
19/26	5	118.00
24/29	5	145.00
29/32	5	168.00



0230 - Receiver Adapters, Plain Bend, Long,
 Made from Heat Resistant Borosilicate Glass, KAVIN

Socket	Qty/pack	Price/Pc
10/19	5	93.00
14/23	5	100.00
19/26	5	118.00
24/29	5	154.00
29/32	5	176.00
34/35	5	193.00



0240 - Receiver Adapters, Straight Delivery,
 Made from Heat Resistant Borosilicate Glass, KAVIN

Socket	Qty/pack	Price/Pc
14/23	5	100.00
19/26	5	118.00
24/29	5	154.00
29/32	5	176.00
34/35	5	200.00
45/40	5	235.00




0250 - Receiver Adapters, Plain Bend, Socket to Cone,
 Made from Heat Resistant Borosilicate Glass, KAVIN

Socket	Cone	Qty/pack	Price/Pc
14/23	14/23	5	160.00
19/26	19/26	5	184.00
24/29	24/29	5	184.00
19/26	24/29	5	245.00
29/32	29/32	5	345.00
34/35	34/35	5	470.00





0260 - Receiver Adapters, Bent with vent,
Made from Heat Resistant Borosilicate Glass, KAVIN



Socket	Cone	Qty/pack	Price/Pc
10/19	10/13	4	143.00
14/23	14/23	4	150.00
19/26	19/26	4	168.00
19/26	24/29	4	168.00
24/29	24/29	4	227.00
29/32	29/32	4	345.00

0270 - Receiver Adapters, Bent with Vacuum Connection,
Made from Heat Resistant Borosilicate Glass, KAVIN



Socket	Cone	Qty/pack	Price/Pc
10/19	10/13	4	168.00
14/23	14/23	4	210.00
19/26	19/26	4	250.00
19/26	24/29	4	268.00
24/29	24/29	4	294.00
29/32	29/32	4	378.00



0280 - Receiver Adapters, Vertical with Vacuum Connection,
Made from Heat Resistant Borosilicate Glass, KAVIN

Socket	Cone	Qty/pack	Price/Pc
19/26	24/29	2	250.00
14/23	24/29	2	268.00



0290 - Receiver Adapters, Bent with Socket Connection,
Made from Heat Resistant Borosilicate Glass, KAVIN

Socket	Cone	Qty/pack	Price/Pc
14/23	14/23	2	218.00
19/26	19/26	2	250.00
19/26	24/29	2	268.00
24/29	24/29	2	300.00



0300 - Receiver Adapters, Bent with Multiple Connection (PIG),
Made from Heat Resistant Borosilicate Glass, KAVIN

Socket	Cone	Qty/pack	Price/Pc
10/19	10/13	2	840.00
14/23	14/23	2	925.00
19/26	14/23	2	925.00
29/32	14/23	2	965.00



0310 - Receiver Adapters, Multi-Limb, Vacuum,
Made from Heat Resistant Borosilicate Glass, KAVIN

Socket	Cone	Qty/pack	Price/Pc
19/26	14/23	2	2925.00



0320 - Receiver Adapters, Vertical Delivery,
Made from Heat Resistant Borosilicate Glass, KAVIN

Socket	Cone	Qty/pack	Price/Pc
19/26	19/26	2	184.00



0330 - Receiver Adapters, Sloping Delivery,
Made from Heat Resistant Borosilicate Glass, KAVIN

Socket	Cone	Qty/pack	Price/Pc
19/26	19/26	2	184.00



0340 - Intermediate Receiver Adapters (PERKIN)
Made from Heat Resistant Borosilicate Glass, KAVIN

Socket	Cone	Cap.ml	Qty/pack	Price/Pc
14/23	14/23	10	1	1840.00
19/26	19/26	32	1	1840.00
19/26	24/29	32	1	2000.00
24/29	24/29	32	1	2180.00



0350 - Recovery Bend, Sloping
Made from Heat Resistant Borosilicate Glass, KAVIN

Flask Cone	Cond. Cone	Qty/pack	Price/Pc
14/23	14/23	4	135.00
19/26	19/26	4	168.00
24/29	14/23	4	210.00
24/29	19/26	4	210.00
24/29	24/29	4	210.00
29/32	19/26	4	210.00
29/32	29/32	4	260.00
34/35	19/26	4	335.00
34/35	24/29	4	420.00
34/35	34/35	4	430.00




0360 - Recovery Bend, Vertical,
Made from Heat Resistant Borosilicate Glass, KAVIN

Flask Cone	Cond. Cone	Qty/pack	Price/Pc
14/23	14/23	2	210.00
19/26	19/26	2	250.00
24/29	19/26	2	294.00
24/29	24/29	2	298.00
34/35	34/35	2	1090.00




0370 - Plain Stillheads, Sloping,
Made from Heat Resistant Borosilicate Glass, KAVIN

Flask Cone	Cond. Cone	Qty/pack	Price/Pc
10/13	10/13	2	184.00
14/23	14/23	2	193.00
19/26	19/26	2	218.00
19/26	19/26	2	235.00
24/29	19/26	2	252.00
24/29	24/29	2	268.00
29/32	19/26	2	268.00
29/32	29/32	2	360.00
34/35	19/26	2	378.00
34/35	24/29	2	386.00
34/35	34/35	2	1090.00





0380 - Claisen Heads, Sloping,
Made from Heat Resistant Borosilicate Glass, KAVIN



Flask Cone	Cond.Cone	Socket	Qty/pack	Price/Pc
14/23	14/23	14/23	2	268.00
19/26	19/26	14/23	2	268.00
24/29	19/26	14/23	2	275.00
24/29	24/29	14/23	2	420.00
29/32	29/32	14/23	2	470.00
34/35	19/26	14/23	2	479.00

0390 - Splash Heads, Adapters,
Made from Heat Resistant Borosilicate Glass, KAVIN




Socket	Cone	Qty/pack	Price/Pc
19/26	24/29	2	325.00
24/29	24/29	2	335.00

0400 - Splash Heads, Sloping,
Made from Heat Resistant Borosilicate Glass, KAVIN




Flask Cone	Cone	Qty/pack	Price/Pc
19/26	19/26	2	310.00
24/29	19/26	2	335.00
24/29	24/29	2	335.00
29/32	19/26	2	345.00
29/32	29/32	2	440.00
34/35	19/26	2	385.00
34/35	24/29	2	395.00
34/35	34/35	2	1090.00

0410 - Splash Heads, Vertical,
Made from Heat Resistant Borosilicate Glass, KAVIN




Flask Cone	Cone	Qty/pack	Price/Pc
19/26	19/26	2	345.00
24/29	19/26	2	378.00
24/29	24/29	2	445.00

0420 - Steam Distillation Heads, Sloping,
Made from Heat Resistant Borosilicate Glass, KAVIN



Flask Cone	Cone	Qty/pack	Price/Pc
24/29	19/26	2	435.00
34/35	19/26	2	453.00

0430 - Splash Head, Double Bulb, Sloping,
Made from Heat Resistant Borosilicate Glass, KAVIN



Flask Cone	Cond.Cone	Socket	Qty/pack	Price/Pc
24/29	19/26	14/23	1	1930.00

0440 - Anti-Foaming, Splash Head, Sloping
Made from Heat Resistant Borosilicate Glass, KAVIN



Flask Cone	Cond.Cone	Socket	Qty/pack	Price/Pc
24/29	19/26	14/23	1	2520.00



0450 - Thermometer Pocket, Closed End,
Made from Heat Resistant Borosilicate Glass, KAVIN

Cone	Stem Length mm	Qty/pack	Price/Pc
14/23	42	2	100.00

0460 - Thermometer Pocket, Open end for Closing
Made from Heat Resistant Borosilicate Glass, KAVIN

Cone	Stem Length mm	Qty/pack	Price/Pc
14/23	200	2	100.00
14/23	400	2	135.00
19/26	200	2	140.00
19/26	400	2	150.00

0470 - Air / Steam Inlet Tubes, Straight,
Made from Heat Resistant Borosilicate Glass, KAVIN

Cone	Stem Length mm	Qty/pack	Price/Pc
7/16	100	2	84.00
10/19	100	2	93.00
14/23	150	2	100.00
19/26	150	2	117.00
24/29	150	2	126.00
29/32	150	2	140.00

0480 - Air / Steam Inlet Tubes, Bend,
Made from Heat Resistant Borosilicate Glass, KAVIN

Cone	Stem Length mm	Qty/pack	Price/Pc
10/19	100	2	93.00
14/23	150	2	100.00
19/26	150	2	117.00
24/29	150	2	126.00

0490 - Drying Tubes, Angled,
Made from Heat Resistant Borosilicate Glass, KAVIN


Cone	Qty/pack	Price/Pc
10/19	2	100.00
14/23	2	117.00
19/26	2	125.00
24/29	2	135.00
29/32	2	160.00


0500 - Tube Adapters for Vacuum Release,
Made from Heat Resistant Borosilicate Glass, KAVIN

Cone	Qty/pack	Price/Pc
19/26	2	525.00
24/29	2	525.00

0510 - Tube Adapters for Vacuum Control,
Made from Heat Resistant Borosilicate Glass, KAVIN

Cone	Qty/pack	Price/Pc
19/26	2	525.00
24/29	2	525.00





02026 - Pycnometers,

Made from Heat Resistant Borosilicate Glass, KAVIN

Capacity ml	Qty/pack	Price/Pc
10	4	450.00
25	4	500.00
50	4	525.00

1001-Beakers, Low form,with Graduation & spout,

Made from Heat Resistant Borosilicate Glass, KAVIN



Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
5	20	25	10	35.00
10	25	32	10	37.00
25	34	45	10	39.00
50	40	55	10	42.00
100	50	73	10	44.00
150	55	80	10	53.00
250	65	98	6	62.00
400	80	110	6	95.00
500	85	120	6	98.00
600	90	125	6	133.00
1000	105	145	6	195.00
2000	130	190	4	450.00

1003-Beakers, Tall form with graduation & spout,

Made from Heat Resistant Borosilicate Glass, KAVIN



Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
50	38	70	12	48.00
100	48	80	12	52.00
150	54	95	12	61.00
250	60	120	6	72.00
400	70	130	6	105.00
500	75	140	6	106.00

1100 - Bottles, B.O.D with I/C Stopper,

Made from Heat Resistant Borosilicate Glass, KAVIN



Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
125	55	120	12	245.00
300	70	150	12	325.00

1110 - Bottles,Reagent Narrow Mouth with I/C Stopper,

Made from Heat Resistant Borosilicate Glass, KAVIN



Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
60	45	95	12	140.00
125	55	95	12	175.00
250	65	150	6	190.00
500	80	170	6	275.00
1000	100	200	2	435.00
2000	130	265	1	675.00

**1120 - Bottles,Reagent Narrow Mouth with I/C Stopper, AMBER**

Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
60	45	95	12	205.00
125	55	95	12	265.00
250	65	150	6	330.00
500	80	170	6	475.00
1000	100	200	2	710.00
2000	130	265	1	900.00

**1130 - Bottles,Reagents Wide Mouth with I/C Stopper,**

Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
60	45	95	12	173.00
125	55	95	12	210.00
250	65	150	6	250.00
500	80	170	6	335.00
1000	100	200	2	485.00
2000	130	265	1	725.00

**1140 - Bottles,Reagent Wide Mouth with I/C Stopper, AMBER**

Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
60	45	95	12	225.00
125	55	95	12	255.00
250	65	150	6	320.00
500	80	170	6	575.00
1000	100	200	2	775.00
2000	130	265	1	820.00

**1150 - Bottles, Laboratory Graduated,With Screwcap and Pouring Ring,**

Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
25	36	74	10	130.00
50	46	91	10	160.00
100	56	105	10	245.00
250	70	143	10	263.00
500	86	181	10	340.00
1000	101	230	10	520.00
2000	136	267	6	760.00

**1160 - Bottles, Laboratory Graduated,With Screwcap and Pouring Ring, AMBER**

Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
25	36	74	10	175.00
50	46	91	10	195.00
100	56	105	10	283.00
250	70	143	10	425.00
500	86	181	10	520.00
1000	101	230	10	650.00
2000	136	267	10	1000.00






1170 - Bottles, Specific Gravity with I/C Capillary bore Glass Stopper,
Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
10	24	45	10	102.00
25	38	55	10	122.00
50	49	65	6	135.00



1180 - Bottles, Specific Gravity with I/C Capillary bore Teflon Stopper,
Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
10	24	45	10	220.00
25	38	55	10	260.00
50	49	65	6	295.00



1190 - Washbottles Complete :
Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Qty/pack	Price/Pc
250	4	240.00
500	4	362.00
1000	2	410.00



1210 - Bottles, Weighing Tall from with I/C Stopper,
Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
5	20	40	10	79.00
15	25	60	10	90.00
25	30	65	10	123.00
60	40	90	6	165.00



1220 - Bottles, Weighing Squat from with I/C Stopper,
Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
20	50	35	10	230.00
40	65	40	10	274.00



1310 - Burettes, Straight bore Glass Solid Stopback : Class B,
Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Grad.Civ	Tol \pm mm	Qty/pack	Price/Pc
10	0.05	0.05	10	180.00
25	0.10	0.10	10	190.00
50	0.10	0.10	10	205.00
100	0.10	0.10	4	220.00



1320 - Burettes, Rotoflow Telfon Screwcock : Class B,
Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Grad.Civ	Tol \pm mm	Qty/pack	Price/Pc
10	0.05	0.05	10	180.00
25	0.10	0.10	10	190.00
50	0.10	0.10	10	205.00
100	0.10	0.10	4	220.00



1330 - Burettes, Straight bore Teflon Stopcock : Class B,
 Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Grad.Div	Tol \pm mm	Qty/pack	Price/Pc
10	0.05	0.05	10	365.00
25	0.10	0.10	10	380.00
50	0.10	0.10	10	397.00
100	0.10	0.10	4	425.00



1340 - Burettes, Straight bore Glass Stopcock : Class B - Schellbach,
 Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Grad.Div	Tol \pm mm	Qty/pack	Price/Pc
10	0.05	0.05	10	198.00
25	0.10	0.10	10	240.00
50	0.10	0.10	10	295.00
100	0.10	0.10	4	320.00



1350 - Burettes, rotoflow Telfon Screw Cock : Class B - Schellbach,
 Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Grad.Div	Tol \pm mm	Qty/pack	Price/Pc
10	0.05	0.05	10	198.00
25	0.10	0.10	10	240.00
50	0.10	0.10	10	295.00
100	0.10	0.10	4	320.00



1360 - Burettes, Straight bore Telfon Stopcock : Class B - Schellbach,
 Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Grad.Div	Tol \pm mm	Qty/pack	Price/Pc
10	0.05	0.05	10	475.00
25	0.10	0.10	10	495.00
50	0.10	0.10	10	545.00
100	0.10	0.10	4	690.00



1370 - Burettes, Automatic, With Reservoir, Straight bore Glasskey Stopcock Class B,
 Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Qty/pack	Price/Pc
10	1	2300.00
25	1	2400.00
50	1	2650.00
100	1	2900.00



2010 - Condensers, Air with I/C inner and outer Joint,
 Made from Heat Resistant Borosilicate Glass, KAVIN

App.Length of Jacket mm	Joint Size	Qty/pack	Price/Pc
100	10/19	4	130.00
200	14/23	4	200.00
200	19/26	4	200.00
300	19/26	4	250.00





2020 - Condensers, Liebig with I/C inner and outer Joint,

Made from Heat Resistant Borosilicate Glass, KAVIN



App.Length of Jacket mm	Joint Size	Qty/pack	Price/Pc
200	B19	1	240.00
300	B24	1	275.00
400	B24	1	345.00

2030 - Condensers, Allihn with I/C inner and outer Joint,

Made from Heat Resistant Borosilicate Glass, KAVIN



App.Length of Jacket mm	Joint Size	Qty/pack	Price/Pc
200	B19	1	295.00
300	B24	1	330.00
400	B24	1	430.00

2040 - Condensers, Graham Coil with I/C inner and outer Joint,

Made from Heat Resistant Borosilicate Glass, KAVIN



App.Length of Jacket mm	Joint Size	Qty/pack	Price/Pc
300	B24	1	460.00
400	B24	1	570.00
500	B24	1	650.00

2050 - Condensers, Double Surface Davies with I/C inner and outer Joint,


Made from Heat Resistant Borosilicate Glass, KAVIN



App.Length of Jacket mm	Joint Size	Qty/pack	Price/Pc
150	B19	1	490.00
200	B19	1	945.00
300	B19	1	1100.00
300	B24	1	1350.00
300	B29	1	1480.00

2060 - Condensers, Friedrichs, Drip Tip, Interchangeable Inner Joint and Interchangeable Outer Joint,

Made from Heat Resistant Borosilicate Glass, KAVIN



Approx Overall Height mm	Interchangeable		Qty pack	Price Pc
	Inner	Outer		
350	24/29	24/29	5	1890.00



3010 - Cylinder Nessler, Colour Comparison

Made from Heat Resistant Borosilicate Glass KAVIN

Cap.ml	Qty/pack	Price/Pc
50ml	10	68.00
100ml	10	95.00



3020 - Cylinder Measuring Class B, Round base

Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Grad.Civ	Tol \pm mm	Qty/pack	Price/Pc
5	0.1	0.1	10	55.00
10	0.2	0.2	10	60.00
25	0.5	0.5	10	85.00
50	1.0	1.0	10	103.00
100	1.0	1.0	10	142.00
250	2.0	2.0	10	270.00
500	5.0	5.0	10	330.00
1000	10.0	10.0	4	550.00
2000	20.0	10.0	2	950.00



3030 - Cylinder Measuring Class B, Hexagonal base

Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Grad.Div	Tol \pm mm	Qty/pack	Price/Pc
5	0.1	0.1	10	130.00
10	0.2	0.2	10	140.00
25	0.5	0.5	10	160.00
50	1.0	1.0	10	180.00
100	1.0	1.0	10	215.00
250	2.0	2.0	10	380.00
500	5.0	5.0	10	510.00
1000	10.0	10.0	4	825.00
2000	20.0	10.0	2	1430.00



3040 - Cylinder Measuring Class B, Hexagonal base, Schellbach

Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Grad.Div	Tol \pm mm	Qty/pack	Price/Pc
5	0.1	0.1	10	153.00
10	0.2	0.2	10	160.00
25	0.5	0.5	10	183.00
50	1.0	1.0	10	225.00
100	1.0	1.0	10	245.00
250	2.0	2.0	10	440.00
500	5.0	5.0	10	590.00
1000	10.0	10.0	4	895.00
2000	20.0	10.0	2	1500.00





3050 - Crow Receivers, Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Tolerance+ml	Qty/pack	Price/Pc
25	0.1 & 0.2	4	500.00
50	0.2 & 0.4	4	545.00
100	0.4 & 0.8	4	580.00

2110 - Colomns, Chromatography, Plain With PTFE Stopcock, Made from Heat Resistant Borosilicate Glass, KAVIN



Capacity ml	Bore mm	Qty pack	Price Pc
200	10	5	200.00
300	10	5	220.00
300	18	5	240.00
450	30	5	310.00
500	10	2	245.00
500	18	2	268.00
600	30	2	400.00
600	40	2	545.00
1000	40	1	1090.00

2120 - Colomns, Chromatography, Plain With Boroflo Screwcock Made from Heat Resistant Borosilicate Glass, KAVIN



Capacity ml	Bore mm	Qty pack	Price Pc
200	10	5	570.00
300	10	5	580.00
300	18	5	655.00
450	30	5	840.00
500	10	2	755.00
500	18	2	925.00
600	30	2	1050.00
600	40	2	1160.00
1000	40	1	1340.00

2130 - Colomns, Chromatography, Plain With Glass Stopcock Made from Heat Resistant Borosilicate Glass, KAVIN



Capacity	Bore	Qty	Price
200	10	5	190.00
300	10	5	190.00
300	18	5	215.00
450	30	5	330.00
500	10	2	259.00
500	18	2	285.00
600	30	2	465.00
600	40	2	580.00
1000	40	1	1160.00



2140 - Colomns, Chromatography, Plain With Sintered Disc and Glass Stopcock,

Made from Heat Resistant Borosilicate Glass, KAVIN


Extractor Length	Bore	Qty	Price
200	10	5	350.00
300	10	5	380.00
300	18	5	480.00
450	30	5	410.00
500	10	2	440.00
500	18	2	515.00
600	30	2	610.00
600	40	2	980.00
1000	40	1	1350.00



2150 - Colomns, Chromatography, Plain With Sintered Disc and PTFE Stopcock,

Made from Heat Resistant Borosilicate Glass, KAVIN

Extractor Length	Bore	Qty	Price
200	10	5	670.00
300	10	5	670.00
300	18	5	980.00
450	30	5	1250.00
500	10	2	890.00
500	18	2	1050.00
600	30	2	1350.00
600	40	2	1695.00
1000	40	1	2100.00



2160 - Colomns, Chromatography, Plain With Sintered Disc and Boroflo Stopcock,

Made from Heat Resistant Borosilicate Glass, KAVIN

Extractor Length	Bore	Qty	Price
200	10	5	310.00
300	10	5	320.00
300	18	5	400.00
450	30	5	500.00
500	10	2	405.00
500	18	2	425.00
600	30	2	560.00
600	40	2	980.00
1000	40	1	1350.00





2170 - Fractionating Column, Vigrex

Made from Heat Resistant Borosilicate Glass, KAVIN



Socket	Cone	Length mm	Qty/pack	Price/Pc
24/29	24/29	360	2	840.00
24/29	24/29	600	2	1500.00
29/32	39/32	360	2	1130.00
29/32	39/32	600	2	1680.00

2190 - Spray for Chromatography

Made from Heat Resistant Borosilicate Glass, KAVIN



Cap.ml	Qty/pack	Price/Pc
25	2	550.00
50	2	550.00
100	2	660.00

3110 - Dessiccators, Plain with Cover Knob Top Kavin, Toughened Natural Glass



D mm	Qty/pack	Price/Pc
150	1	2300.00
210	1	2550.00
240	1	3300.00
300	1	5100.00

3120 - Dessiccators, Vaccum with Cover Knob Top Kavin, Toughened Natural Glass



D mm	Qty/pack	Price/Pc
150	1	4200.00
210	1	4400.00
240	1	5800.00
300	1	8500.00

3130 - Dishes Petri, Made from Heat Resistant Borosilicate Glass, KAVIN



Approx O.D.xht mm	Qty/pack	Price/Pc
50x17	20	80.00
80x17	10	95.00
100x17	10	105.00
150x20	10	320.00

3140 - Dishes, Crystallizing,

Made from Heat Resistant Borosilicate Glass, KAVIN



Approx O.D.xht mm	Qty/pack	Price/Pc
100x50	10	265.00
150x75	10	460.00
190x100	4	770.00

3150 - Dishes, Evaporating, Flat Bottom, With Pour Out,

Made from Heat Resistant Borosilicate Glass, KAVIN



Approx O.D.xht mm	Qty/pack	Price/Pc
80x45	60	290.00
105x55	50	460.00
150x80	20	750.00
200x100	8	1600.00



3210 - Dean & Stark Apparatus,

Made from Heat Resistant Borosilicate Glass, KAVIN

Specefication	Qty/pack	Price/Pc
Receiver 10ml, R.B. Flask 1000 ml. Liebig Condenser,	1	1230.00



3220 - Extraction Apparatus, Soxhlet Complete with allihn Condenser and F.B. Flasks

Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Flask Cap.ml	Qty/pack	Price/Pc
60	150	1	800.00
100	250	1	1200.00
200	500	1	1625.00



3230 - Extraction Apparatus, Spare Allihn Condenser

Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Joint Bottom	Qty/pack	Price/Pc
For 60	34/35	1	315.00
For 100	40/38	1	430.00
For 200	50/42	1	600.00



3240 - Extraction Apparatus, Spare Extractor

Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	I/C Joint Top	I/C Joint Bottom	Qty/pack	Price/Pc
60	34/35	24/29	1	370.00
100	40/38	24/29	1	550.00
200	50/42	24/29	1	700.00



4010 - Flasks Conical, Narrow Mouth

Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
10	34	60	20	45.00
25	42	70	20	46.00
50	51	82	20	47.00
100	64	105	10	55.00
150	74	120	10	70.00
250	84	140	6	83.00
500	104	175	6	120.00
1000	130	220	6	195.00
2000	163	285	1	460.00



4020 - Flasks Conical, wide Mouth with Graduation

Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
100	64	105	10	63.00
250	72	124	6	87.00
500	104	175	6	135.00
1000	131	210	6	210.00





4030 - Flasks Conical, with I/C Stopper

Made from Heat Resistant Borosilicate Glass, KAVIN



Cap.ml	I/C Stopper	Qty/pack	Price/Pc
100	B19	10	125.00
100	B24	10	130.00
250	B24	6	140.00
250	B29	6	145.00
500	B24	4	225.00
500	B29	4	270.00
1000	B24	2	350.00
1000	B29	2	385.00
2000	B24	1	650.00
2000	B29	1	695.00

4210 - Flasks Conical, With Screw Cap,

Made from Heat Resistant Borosilicate Glass, KAVIN



Cap.ml	Approx O.D.x Height mm	Qty/pack	Price/Pc
100	64x110	10	195.00
150	72x125	10	221.00
250	85x145	10	273.00
500	104x180	10	336.00
1000	131x225	10	641.00

4040 - Flasks Iodine,with I/C Stopper

Made from Heat Resistant Borosilicate Glass, KAVIN



Cap.ml	I/C Stopper	Qty/pack	Price/Pc
100	B24	10	125.00
250	B24	6	170.00
500	B24	4	235.00
1000	B24	2	380.00

4050 - Flasks Flat Bottom, Narrow Neck Plain

Made from Heat Resistant Borosilicate Glass, KAVIN



Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
25	42	90	10	48.00
50	51	95	10	50.00
100	64	110	10	60.00
150	75	137	10	74.00
250	85	145	10	75.00
500	105	178	6	115.00
1000	133	212	6	170.00
2000	167	230	2	333.00

4060 - Flasks Round Bottom, Narrow Neck Plain

Made from Heat Resistant Borosilicate Glass, KAVIN



Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
25	42	90	10	48.00
50	51	95	10	50.00
100	64	110	10	60.00
150	75	137	10	74.00
250	85	145	10	75.00
500	105	178	6	115.00
1000	133	212	6	170.00
2000	167	230	2	333.00



4070 - Flasks, Flat Bottom, Centre Neck with I/C Joint
 Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Joint Size	Qty/pack	Price/Pc
5	14/23	10	50.00
10	14/23	10	70.00
25	14/23	10	80.00
25	19/ 23	10	80.00
50	19/26	10	100.00
50	24/29	10	100.00
50	29/32	10	100.00
50	34/35	10	100.00
100	19/26	10	120.00
100	24/29	10	120.00
100	29/32	10	120.00
100	34/35	10	120.00
150	19/26	8	135.00
150	24/29	8	135.00
150	29/32	8	135.00
150	34/35	8	135.00
250	19/26	8	175.00
250	24/29	8	175.00
250	29/32	8	175.00
250	34/35	8	175.00
500	19/26	4	225.00
500	24/29	4	225.00
500	29/32	4	225.00
500	34/35	4	225.00
1000	19/26	4	360.00
1000	24/29	4	360.00
1000	29/32	4	375.00
1000	34/35	4	392.00



4080 - Flasks, Round Bottom, Centre Neck with I/C Joint
 Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Joint Size	Qty/pack	Price/Pc
5	14/23	10	50.00
10	14/23	10	70.00
25	14/23	10	80.00
25	19/23	10	80.00
50	19/26	10	100.00
50	24/29	10	100.00
50	29/32	10	100.00
50	34/35	10	100.00
100	19/26	10	120.00
100	24/29	10	120.00
100	29/32	10	120.00





Cap.ml	Joint Size	Qty/pack	Price/Pc
100	34/35	10	120.00
150	19/26	8	135.00
150	24/29	8	135.00
150	29/32	8	135.00
150	34/35	8	135.00
250	19/26	8	175.00
250	24/29	8	175.00
250	29/32	8	175.00
250	34/35	8	175.00
500	19/26	4	225.00
500	24/29	4	225.00
500	29/32	4	225.00
500	34/35	4	225.00
1000	19/26	4	290.00
1000	24/29	4	290.00
1000	29/32	4	290.00
1000	34/35	4	290.00
2000	19/26	2	550.00
2000	24/29	2	550.00
2000	29/32	2	550.00
2000	34/35	2	550.00

4090 - Flasks, Round Bottom, Two Necks, Centre Neck and One Side Angled Neck with I/C Joint. Made from Heat Resistant Borosilicate Glass, KAVIN



Cap.ml	Centre Neck	Side Neck	Qty/pack	Price/Pc
100	24/29	14/23	5	175.00
100	24/29	19/26	5	175.00
250	24/29	19/26	4	220.00
250	29/32	14/23	4	220.00
250	34/35	19/26	4	220.00
500	19/26	14/23	2	250.00
500	24/29	19/26	2	250.00
1000	29/32	14/23	1	325.00
1000	34/35	19/26	1	325.00

4110 - Flasks, Round Bottom, Two Necks, Centre Neck and Parallel Neck with I/C Joint. Made from Heat Resistant Borosilicate Glass, KAVIN



Cap.ml	Centre Neck	Side Neck	Qty/pack	Price/Pc
100	24/29	14/23	2	260.00
100	24/29	19/26	2	260.00
250	29/32	14/23	2	360.00
250	34/35	19/26	2	370.00
500	19/26	14/23	2	430.00
500	24/29	19/26	2	450.00
1000	29/32	14/23	1	325.00
1000	34/35	19/26	1	325.00
2000	29/32	14/23	1	650.00
2000	34/35	19/26	1	650.00



4120 - Flasks, Round Bottom, Three Necks, Centre Neck and Two Angled Neck with I/C Joint
 Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Centre Neck	Side Neck	Qty/pack	Price/Pc
100	14/23	14/23	2	229.00
100	19/26	14/23	2	240.00
100	24/29	14/23	2	240.00
100	24/29	19/26	2	245.00
250	19/26	14/23	2	290.00
250	19/26	19/26	2	290.00
250	24/29	14/23	2	290.00
250	24/29	19/26	2	290.00
250	29/32	19/26	2	375.00
500	19/26	14/23	2	375.00
500	19/26	19/26	2	375.00
500	24/29	14/23	2	375.00
500	24/29	19/26	2	375.00
500	29/32	14/23	2	390.00
500	29/32	19/26	2	390.00
1000	24/29	14/23	1	610.00
1000	24/29	19/26	1	610.00
1000	34/35	19/26	1	640.00
2000	24/29	14/23	1	764.00
2000	24/29	19/26	1	764.00



4130- Flasks, Filter Conical Shape with Glass Hose Connection

Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Approx O.D	Ht mm	Qty/pack	Price/Pc
100	64	105	10	131.00
250	85	150	6	240.00
500	105	184	4	310.00
1000	136	245	2	570.00
2000	168	295	1	1350.00



4190 - Flasks Culture, Haffkine,

Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Approx O.D	Qty/pack	Price/Pc
3000	245x273	2	2050.00
4000	246x310	2	2400.00



4140 - Flasks, Distilling with side arm

Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Approx O.D	Ht mm	Qty/pack	Price/Pc
100	64	210	6	140.00
250	68	214	4	200.00
500	88	226	2	225.00
1000	110	250	1	340.00



4150 - Flasks, Kjeldahi Plain

Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Approx O.D	Ht mm	Qty/pack	Price/Pc
300	97	335	6	200.00
500			6	245.00
800	110	375	6	325.00





4160 - Flasks, Kjeldahi with i/C Joint,

Made from Heat Resistant Borosilicate Glass, KAVIN



Cap.ml	Joint Size	Qty/pack	Price/Pc
300	B24	6	300.00
500	B24	6	335.00
800	B24	6	435.00

4170 - Flasks, Volumetric with I/C Stopper Class B

Made from Heat Resistant Borosilicate Glass, KAVIN



Cap.ml	Tel±mm	Joint	Qty/pack	Price/Pc
5	0.05	B10	10	75.00
10	0.05	B10	10	80.00
20	0.08	B10	10	85.00
25	0.08	B10	10	92.00
50	0.12	B10	10	99.00
100	0.20	B14	10	125.00
200	0.30	B14	10	146.00
250	0.30	B14	6	155.00
500	0.50	B19	6	210.00
1000	0.80	B19	6	315.00

4180 - Flasks, Volumetric with I/C Stopper Class B

Made from Heat Resistant Borosilicate Glass, KAVIN


AMBER



Cap.ml	Tel±mm	Joint	Qty/pack	Price/Pc
5	0.05	B10	10	103.00
10	0.05	B10	10	123.00
20	0.08	B10	10	129.00
25	0.08	B10	10	133.00
50	0.12	B10	10	145.00
100	0.20	B14	10	215.00
200	0.30	B14	10	225.00
250	0.30	B14	6	245.00
500	0.50	B19	6	378.00
1000	0.80	B19	6	593.00

4190 - Flasks, Boiling, Pear Shaped, Short Neck, With Interchangeable Joint

Made from Heat Resistant Borosilicate Glass, KAVIN



Capacity ml	Interchangeable	Qty/pack	Price/Pc
5	10/19	10	80.00
5	14/23	10	80.00
10	10/19	10	95.00
10	14/23	10	95.00
25	14/23	5	140.00
50	14/23	10	160.00
100	24/29	10	235.00



4210 - Flasks, Pear Shape Suitable for Rotary Evaporators with Ground Socket,
 Made from Heat Resistant Borosilicate Glass, KAVIN

Capacity ml	Ground Joint	Qty/pack	Price/Pc
50	24/29	5	300.00
50	29/32	5	300.00
100	24/29	5	390.00
100	29/32	5	390.00
250	24/29	5	450.00
250	29/32	5	450.00
500	24/29	5	630.00
500	29/32	5	630.00
1000	24/29	5	940.00
1000	29/32	5	940.00
2000	24/29	5	1400.00
2000	29/32	5	1400.00



4310 - Funnels, 60 Long Stem,
 Made from Heat Resistant Borosilicate Glass, KAVIN

Dia. mm	Qty/pack	Price/Pc
50	6	42.00
65	6	53.00
75	6	67.00
100	2	97.00



4320 - Funnels, 60 Short Stem, for Powder,
 Made from Heat Resistant Borosilicate Glass, KAVIN

Dia. mm	Qty/pack	Price/Pc
125	4	159.00
150	4	199.00
200	2	388.00



4330 - Funnels, Separating Pear Shape with Rotaflow Teflon Screwcock & I/C Stopper
 Made from Heat Resistant Borosilicate Glass, KAVIN

Dia. mm	Qty/pack	Price/Pc
60	4	160.00
125	4	205.00
250	4	249.00
500	2	380.00
1000	2	792.00



4340 - Funnels, Separating Pear Shape with Glass Stopcock & I/C Stopper
 Made from Heat Resistant Borosilicate Glass, KAVIN

Dia. mm	Qty/pack	Price/Pc
60	4	130.00
125	4	175.00
250	4	225.00
500	2	310.00
1000	2	475.00



2140 - Funnels, Separating Pear Shape with Teflon Stopcock & I/C Stopper
 Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Qty/pack	Price/Pc
60	4	495.00
125	4	523.00
250	4	555.00
500	2	745.00
1000	2	1100.00





4350 - Funnels, Separating Globe Shape with Rotoflow Teflon Screwcock & I/C Stopper
 Made from Heat Resistant Borosilicate Glass, KAVIN



Dia. mm	Qty/pack	Price/Pc
60	4	495.00
125	4	523.00
250	4	555.00
500	2	745.00

4360 - Funnels, Separating Globe Shape with Glass Stopcock & I/C Stopper
 Made from Heat Resistant Borosilicate Glass, KAVIN



Dia. mm	Qty/pack	Price/Pc
60	4	145.00
125	4	195.00
250	4	235.00
500	2	365.00
1000	2	725.00

4370 - Funnels, Separating Globe Shape with Teflon Stopcock & I/C Stopper
 Made from Heat Resistant Borosilicate Glass, KAVIN



Dia. mm	Qty/pack	Price/Pc
60	4	495.00
125	4	523.00
250	4	525.00
500	2	745.00
1000	2	1100.00

4390 - Pressure Equalising Funnels, Spherical,
 Made from Heat Resistant Borosilicate Glass, KAVIN



Capacity ml	Socket	Cone	Qty/pack	Price/Pc
100	19/26	19/26	2	630.00
100	24/29	19/29	2	630.00
100	24/29	24/29	2	630.00
100	29/32	29/32	2	630.00
250	24/29	19/26	2	850.00
250	24/29	24/29	2	850.00
250	29/32	29/32	2	850.00
250	34/35	34/35	2	850.00
500	34/35	34/35	2	1400.00

4410 - Pressure Equalising Funnels, Cylindrical,
 Made from Heat Resistant Borosilicate Glass, KAVIN



Capacity ml	Socket	Cone	Qty/pack	Price/Pc
100	24/29	19/26	1	945.00
250	24/29	24/29	1	1100.00
500	29/32	29/32	1	1260.00



4420 - Filter Crucibles with Sintered Disc,
 Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Porosity Grade	Qty/pack	Price/Pc
15	G1	10	95.00
15	G2	10	95.00
15	G3	10	95.00
15	G4	10	95.00
30	G1	10	130.00
30	G2	10	130.00
30	G3	10	130.00
30	G4	10	130.00
50	G1	10	175.00
50	G2	10	175.00
50	G3	10	175.00



4430 - Funnels Buchner with Sintered Disc,
 Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Porosity Grade	Qty/pack	Price/Pc
35	G1	2	280.00
35	G2	2	280.00
35	G3	2	280.00
35	G4	2	280.00
80	G1	2	480.00
80	G2	2	480.00
80	G3	2	480.00
80	G4	2	480.00
200	G1	2	720.00
200	G2	2	720.00
200	G3	2	720.00



4110 - Pipettes, Graduated, (Serological Type) Graduated to Tip Class B,
 Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Grad. Interval mm	Tol+ml	Qty/pack	Price/Pc
1	0.01	0.01	10	45.00
1	0.10	0.01	10	45.00
2	0.02	0.02	10	48.00
2	0.10	0.02	10	48.00
5	0.05	0.05	10	55.00
5	0.10	0.10	10	55.00
10	0.10	0.10	10	73.00
25	0.10	0.20	10	120.00



4120 - Pipettes, Graduated, (Mohrs Type) Graduated to Tip Class B,
 Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Grad. Interval mm	Tol+ml	Qty/pack	Price/Pc
1	0.01	0.01	10	45.00
1	0.10	0.01	10	45.00
2	0.02	0.02	10	48.00
2	0.10	0.02	10	48.00
5	0.05	0.05	10	55.00
5	0.10	0.10	10	55.00
10	0.10	0.10	10	73.00
25	0.10	0.20	10	120.00





5120 - Pipettes, Volumetric, Class B,
Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Tol+ml	Qty/pack	Price/Pc
1	0.015	10	45.00
2	0.020	10	45.00
5	0.030	10	48.00
10	0.040	10	48.00
20	0.060	10	55.00
25	0.060	10	55.00
50	0.100	10	73.00
100	0.160	10	120.00

5130 - Pipettes, Artificial Insemination, Cattle, Taper Tip,
Made from Heat Resistant Borosilicate Glass, KAVIN

Length mm	Qty/pack	Price/Pc
450	1	45.00

5140 - Drum Sampling, Pipette,
Made from Heat Resistant Borosilicate Glass, KAVIN

Length mm	Qty/pack	Price/Pc
1000	4	215.00

5150 - Pipettes, Bacteriological, Graduated
Made from Heat Resistant Borosilicate Glass, KAVIN

Capital ml	Tolerance+ml	Qty/pack	Price/Pc
1.1	0.025	10	98.00
2.2	0.04	10	107.00

5160 - Pipettes, Gerber. Milk,
Made from Heat Resistant Borosilicate Glass, KAVIN

Capital ml	Tolerance+ml	Qty/pack	Price/Pc
10.75	0.03	10	98.00

5210 - Tilt Measure, with Ground Reservoir,
Made from Heat Resistant Borosilicate Glass, KAVIN

Capital ml	Reservoir Capacity ml	Qty/pack	Price/Pc
5	250	2	630.00
10	250	2	630.00
20	500	2	700.00
25	500	2	700.00
50	1000	2	780.00

5220 - Tubes, Filter for Ground Crucibles
Made from Heat Resistant Borosilicate Glass, KAVIN

Approx. O.D. x Length mm	Qty/pack	Price/Pc
27 x 150	5	85.00
36 x 160	5	115.00
56 x 170	10	180.00



**6110 - Culture Tubes, Round Bottom with PP Screw Cap,**

Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Qty/pack	Price/Pc
5	100	26.00
15	100	36.00
30	100	42.00
60	100	52.00

**6120 - Culture Tubes, Flat Bottom with PP Screw Cap,**

Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Qty/pack	Price/Pc
5	100	26.00
15	100	36.00
30	100	40.00

**6130 - Culture Tubes, Round Bottom with PP Screw Cap, AMBER**

Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Qty/pack	Price/Pc
5	100	42.00
15	100	72.00
30	100	82.00
60	100	92.00

**6140 - Culture Tubes, Flat Bottom with PP Screw Cap, AMBER**

Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Qty/pack	Price/Pc
5	100	42.00
15	100	72.00
30	100	80.00

**6150 - Centrifuge Tube : Conical Plain**

Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Qty/pack	Price/Pc
15	100	28.00
50	100	62.00

**6160 - Centrifuge Tube : Conical Graduated**

Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Qty/pack	Price/Pc
15	100	34.00
50	100	90.00

**6170 - Blood Sugar Tubes**

Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Graduation	Price/Pc
25	12.5 & 25	44.00
25	4, 12.5 & 25	52.00





6180 - Test Tubes with Rim

Made from Heat Resistant Borosilicate Glass, KAVIN



Size	Qty/pack	Price/Pc
10x75	100	5.00
12x75	100	6.00
12x100	100	8.00
15x125	100	9.00
15x150	100	10.00
18x150	100	11.00
25x150	100	22.00
32x200	25	45.00

6190 - Test Tubes without Rim

Made from Heat Resistant Borosilicate Glass, KAVIN



Size	Qty/pack	Price/Pc
10x75mm	100	6.00
12x75mm	100	7.00
12x100mm	100	8.00
15x125mm	100	9.00
15x150mm	100	10.00
18x150mm	100	11.00
25x150mm	100	20.00
32x200mm	25	42.00

6210 - Solid Stoppers :

Made from Heat Resistant Borosilicate Glass, KAVIN



Size	Qty/pack	Price/Pc
B12	25	35.00
B14	25	40.00
B19	25	45.00
B24	25	75.00
B29	25	95.00
B34	25	135.00
B40	25	190.00
B50	25	200.00
B55	25	210.00

6310 - Membrane filter Holder all glass Assembly

Made from Heat Resistant Borosilicate Glass, KAVIN



Funnel for 47mm dia filter base with sintered disc & Buchner Flask Cap. 1000ml, Funnel 300ml with suitable clamp.	3600.00
---	---------

7010 - Reflux Assembly

Made from Heat Resistant Borosilicate Glass, KAVIN



Components	Qty/pack	Price/Pc
R.B. Flask 100 ml. Liebig Condenser, 300 mm	1	500.00
R.B. Flask 250 ml. Liebig Condenser, 400 mm	1	550.00
R.B. Flask 500 ml. Liebig Condenser, 500 mm	1	630.00



7020 - Reflux Assembly

Made from Heat Resistant Borosilicate Glass, KAVIN

Components	Qty/pack	Price/Pc
R.B. Flask 1000 ml. (24/29 neck) Reversible Condenser,	1	700.00
R.B. Flask 1000 ml. (34/35 neck) Reversible Condenser,	1	780.00



7030 - Reflux Assembly

Made from Heat Resistant Borosilicate Glass, KAVIN


Components	Qty/pack	Price/Pc
R.B. Flask 500 ml. Multiple adapter two necks, stirrer guide Rubber sleeve gland Link stirrer, 450 mm, Vertical Allihn condenser, 300 mm	1	1650.00



7040 - Reaction Assembly,

Made from Heat Resistant Borosilicate Glass, KAVIN


Components	Qty/pack	Price/Pc
Round bottom flask 500 ml, multiple adapter two necks, Stirrer guide, Stuffing Box gland with pulley, Link stirrer 450 mm, Stillhead Thermometer Pocket, Liebig condenser 400 mm, Receiver adapter bent with vent, Round bottom flask 100 ml	1	4000.00



7050 - Reaction Assembly,

Made from Heat Resistant Borosilicate Glass, KAVIN


Components	Qty/pack	Price/Pc
Round bottom flask 1000 ml, multiple adapter two necks, Stirrer guide, Rubber sleeve gland collapsible paddle stirrer 450 mm, Stillhead jacketed coil condenser 150 mm Receiver adapter bent with socket connection, Round bottom flask 250 ml Double surface condenser 200 mm, Cone adapter bent Thermometer with P.T.F.E. Cone	1	4000.00










7060- Reaction Assembly,

Made from Heat Resistant Borosilicate Glass, KAVIN

Components	Qty/pack	Price/Pc
Round bottom flask 500 ml, multiple adapter two necks, Claisen head, Liebig condenser 200 mm, Receiver adapter bent with vent, R.B. flask, 250 ml, Separating funnel, 100 ml, Separating funnel, 100 ml, Stopper, stopper, Thermometer pocket	1	4000.00





	7070- Reaction Assembly, Made from Heat Resistant Borosilicate Glass, KAVIN		
	Components	Qty/pack	Price/Pc
	Round bottom flask 250 ml, multiple adapter two necks, Stillhead Liebig condenser 200 mm, Receiver adapter Thermometer pocket, Seperating funnel, 100 ml, Stopper	1	2350.00
	7080- Distillation Assembly, Made from Heat Resistant Borosilicate Glass, KAVIN		
	Components	Qty/pack	Price/Pc
	Round bottom flask 2000 ml, Splash head, Double surface condenser, 200 mm Receiver adapter	1	2150.00
	7090 - Distillation Assembly, Made from Heat Resistant Borosilicate Glass, KAVIN		
	Components	Qty/pack	Price/Pc
	Round bottom flask 500 ml, Stillhead Liebig condenser 400 mm Receiver adapter, Seperating funnel, 100 ml, Stopper	1	1650.00
			
	7110 - Distillation Assembly, Made from Heat Resistant Borosilicate Glass, KAVIN		
	Components	Qty/pack	Price/Pc
	Round bottom flask 1000 ml, Recovery Bend, Thorpe condensor, 200 mm Receiver adapter	1	2350.00
	7120- Distillation Assembly, Made from Heat Resistant Borosilicate Glass, KAVIN		
	Components	Qty/pack	Price/Pc
	Round bottom flask 250 ml, Recovery Bend, Liebig condensor, 200 mm Receiver adapter	1	850.00
	7130 - Vacuum Distillation Assembly, Made from Heat Resistant Borosilicate Glass, KAVIN		
	Components	Qty/pack	Price/Pc
	Round bottom flask 3000 ml, Multiple adapter, Two necks, Stillhead Jacketed coil condencer, 150 mm, Recovery adapter, bent with socket connection, R.B. Flask, 250ml, Vacuum release stopcock, Thermometer with P.T.F.E. cone, Seperating funnel, 250ml, Stopper	1	5400.00
	7140 - Vacuum Distillation Assembly, Made from Heat Resistant Borosilicate Glass, KAVIN		
	Components	Qty/pack	Price/Pc
	Round bottom flask 100 ml, Multiple adapter, Two necks, Stillhead Air Condenser, 400 mm, Recovery adapter, bent with vent, R.B. Flask, 50 ml Thermometer with P.T.F.E. cone, Seperating funnel, 50ml, Stopper	1	2550.00




**7150 - Vacuum Distillation Assembly,
Made from Heat Resistant Borosilicate Glass, KAVIN**

Components	Qty/pack	Price/Pc
Round bottom flask 500 ml, Splash head, Double bulb, Double Surface condenser, 200 mm Receiver adapter, with Vacuum connection, Round bottom flask, 100 ml Air leak tube	1	2990.00



**7160 - Vacuum Distillation Assembly,
Made from Heat Resistant Borosilicate Glass, KAVIN**

Components	Qty/pack	Price/Pc
Round bottom flask 500 ml, Claisen head, Jacketed coil condenser, 150 mm, Receiver adapter, with Vacuum connection, Round bottom flask, 100 ml, Test tube, 200 ml gas wash bottle fitting, Thermometer with P.T.F.E. Cone, Air leak tube	1	3550.00



**7170 - Vacuum Distillation Assembly,
Made from Heat Resistant Borosilicate Glass, KAVIN**

Components	Qty/pack	Price/Pc
Round bottom flask 1000 ml, Recovery Bend, Double coil condenser, Round Bottom Falsk, Two necks 500 ml Vacuum release stopcock	1	3800.00



**7180 - Vacuum Distillation Assembly,
Made from Heat Resistant Borosilicate Glass, KAVIN**

Components	Qty/pack	Price/Pc
Round bottom flask 1000 ml, Recovery Bend, Jacketed coil condenser, 150 mm, Receiver adapter, bentwith vent, R.B. flask, 250 ml	1	2520.00



**7190 - Kjeldhal Distillation Assembly,
Made from Heat Resistant Borosilicate Glass, KAVIN**

Components	Qty/pack	Price/Pc
Round bottom flask 1000 ml, Multiple adapter ,Two necks, Splash head, vertical, Liebig condenser, 400 mm, Receiver adapter, Seperating funnel, 100 ml	1	2350.00



**7986 - Kjeldhal Distillation Assembly,
Made from Heat Resistant Borosilicate Glass, KAVIN**

Components	Qty/pack	Price/Pc
Kjeldhal flask 800 ml, Multiple adapter ,Two necks, Splash head, vertical, Liebig condenser, 400 mm, Receiver adapter, Seperating funnel, 100 ml	1	2600.00





7220 - Fractionation Assembly,

Made from Heat Resistant Borosilicate Glass, KAVIN



Components	Qty/pack	Price/Pc
Round bottom flask 500 ml, Cup Joint adapter, Ball Joint Adapter, Joint clip, S 35 Air condenser, 300 mm, Multiple adapter, Two necks, Immersion condenser, Stillhead, Thermometer pocket, Receiver adapter, Plain bend Liebig Condenser, 400 mm, Receiver adapter	1	4200.00

7230 - Fractionation Assembly,

Made from Heat Resistant Borosilicate Glass, KAVIN



Components	Qty/pack	Price/Pc
Round bottom flask 250 ml, Air condenser, 400 mm, Stillhead, Thermometer with P.T.F.E. cone, Liebig condenser, 200 mm Bend with multiple connection, 3 x R.B. flask, 25 ml	1	3150.00

7240 - Steam Distillation Assembly,

Made from Heat Resistant Borosilicate Glass, KAVIN



Components	Qty/pack	Price/Pc
R.B. flask, Two necks 500 ml, Splash head, Liebig condenser, 200 mm, Receiver adapter, Air leak tube	1	1650.00

7250 - Steam Distillation Assembly,

Made from Heat Resistant Borosilicate Glass, KAVIN



Components	Qty/pack	Price/Pc
R.B. flask, 500 ml, Stillhead, Air leak tube, Liebig condenser Splash 400 mm, Receiver adapter	1	1600.00

7260 - Steam Distillation Assembly,

Made from Heat Resistant Borosilicate Glass, KAVIN



Components	Qty/pack	Price/Pc
R.B. flask, 500 ml, Steam Distillation head, Liebig condenser 400 mm, Receiver adapter	1	1600.00

7270 - Steam Distillation Assembly,

Made from Heat Resistant Borosilicate Glass, KAVIN



Components	Qty/pack	Price/Pc
R.B. flask, 1000 ml, Multiple adapter, Two necks, Air leak tube, Bent, Recovery bend, Jacketed coil condenser, 150 mm Receiver adapter	1	2350.00

7280 - Recovery Assembly,

Made from Heat Resistant Borosilicate Glass, KAVIN




Components	Qty/pack	Price/Pc
R.B. flask, 250 ml, Multiple adapter, Two necks, Recovery bend, Liebig condenser, 400 mm Receiver adapter Separating funnel, 100 ml	1	2000.00



8010- Recovery Assembly,

Made from Heat Resistant Borosilicate Glass, KAVIN

Components	Qty/pack	Price/Pc
R.B. flask, 2000 ml, Stillhead, Jacketed coil condenser, 150 mm Receiver adapter, Sepaerating Funnel, 100 ml, Stopper	1	3150.00



8020 - Solvent Recovery Assembly,

Made from Heat Resistant Borosilicate Glass, KAVIN

Components	Qty/pack	Price/Pc
R.B. flask, 150ml, Recovery bend, Jacketed coil condenser, 150 mm Receiver adapter	1	2150.00



8030 - Solvent Recovery Assembly,

Made from Heat Resistant Borosilicate Glass, KAVIN


Components	Qty/pack	Price/Pc
R.B. flask, 150ml, Recovery bend, Jacketed coil condenser, 150 mm Receiver adapter	1	1750.00



8040 - Solvent Recovery Assembly,

Made from Heat Resistant Borosilicate Glass, KAVIN


Components	Qty/pack	Price/Pc
R.B. flask, 500ml, Recovery bend, Air Condenser, 200 mm Ether condenser, Buchner flask 250 ml	1	5100.00



8050 - Solvent Recovery Assembly,

Made from Heat Resistant Borosilicate Glass, KAVIN

Components	Qty/pack	Price/Pc
R.B. flask, 500ml, Hopkin Condenser, Receiver adapter	1	1750.00



8060- Solvent Recovery Assembly,

Made from Heat Resistant Borosilicate Glass, KAVIN

Components	Qty/pack	Price/Pc
R.B. flask, 250ml, Plain body extractor, Coil condenser, 150 mm	1	1550.00



8070 - Arsenic Determiration Apparatus,

Made from Heat Resistant Borosilicate Glass, KAVIN

Components	Qty/pack	Price/Pc
Comprising of 100 ml, Bottle and arsenic tube in two parts with rubber cork ends	1	450.00



8080 - Arsenic Determiration Apparatus Gutzeit,

Made from Heat Resistant Borosilicate Glass, KAVIN

Components	Qty/pack	Price/Pc
Comprising of 100 ml, Bottle and inner fitting. The arsenic tube is in two parts with ground flange ends and can be connected together by means of a spring of rubber band	1	1400.00





8110 - Essential Oil Determination Apparatus,
Made from Heat Resistant Borosilicate Glass, KAVIN

Discription	Qty/pack	Price/Pc
Cocking and Middletown - for the determination of volatile oil in drugs-comprising of 1000 ml flask and vertical tube combined with condenser and measuring tube with stopcock. A return tube for the aqueous part of the distillate connects the bottom of the measuring tube and vertical tube.	1	3500.00




8120 - Essential Oil Determination Apparatus,
Made from Heat Resistant Borosilicate Glass, KAVIN

Discription	Qty/pack	Price/Pc
Clevenger-For the determination of volatile oil lighter than water comprising of 1000 ml flask, Oil separatory tube and condenser	1	2350.00



8130 - Essential Oil Determination Apparatus,
Made from Heat Resistant Borosilicate Glass, KAVIN


Discription	Qty/pack	Price/Pc
Clevenger-For the determination of volatile oil lighter than water comprising of 1000 ml flask, Oil separatory tube and condenser	1	2350.00



8140- Distilling Apparatus, with graham Condenser Inter- changeable joint, interchangeable Stopper,

Made from Heat Resistant Borosilicate Glass, KAVIN

Capacity ml	Condenser Jacket mm	Interchangeable	Qty/pack	Price/Pc
500	200	19/26	2	1550.00



8150 - Distilling Apparatus, Ammonia, with graham Condenser Interchangeable joints,
Made from Heat Resistant Borosilicate Glass, KAVIN

Capacity ml	Qty/pack	Price/Pc
500	2	1550.00



8160 - Distilling Apparatus, with Friedrichs Condenser Interchangeable joints, interchangeable Stopper,

Made from Heat Resistant Borosilicate Glass, KAVIN

Capacity ml	Overall Height mm	Interchangeable Joint	Qty/pack	Price/Pc
2000	350	24/29	2	1950.00



8170 - Distilling Apparatus, Compact with Friedrichs Condenser Interchangeable joints,
Made from Heat Resistant Borosilicate Glass, KAVIN

Capacity ml	Overall Height mm	Interchangeable Joint	Qty/pack	Price/Pc
1000	495	34/35	2	1950.00



LIQUID HANDLING PRODUCTS

MICROPIPETTES, Fixed Volume Single channel, Fully Autoclavable, Biopette

Made of High grade engineering plastics, Colour Coded for easy identification, Calibration Conforming to Din 12650 & EN-ISO 8655, Standard ensuring high accuracy & precision, Does not require Calibration after Autoclaving.

Cat No.	Volume μ l	Qty/pack	Price/Pc
F - 2.5	2.5	1	1400.00
F - 5	5	1	1400.00
F - 10	10	1	1400.00
F - 20	20	1	1400.00
F - 25	25	1	1400.00
F - 50	50	1	1400.00
F - 100	100	1	1400.00
F - 200	200	1	1400.00
F - 250	250	1	1400.00
F - 500	500	1	1400.00
F - 1000	1000	1	1400.00
F - 2000	2000	1	1400.00
F - 5000	5000	1	1400.00

MICROPIPETTES, Variable Volume Single channel, Fully Autoclavable, Biopette

Made of High grade engineering plastics, Colour Coded for easy identification, Calibration Conforming to Din 12650 & EN-ISO 8655, Standard ensuring high accuracy & precision, Does not require Calibration after Autoclaving.

Cat No.	Volume μ l	Qty/pack	Price/Pc
V - 2.5	0.1 - 2.5	1	2500.00
V - 10	0.5 - 10	1	2500.00
V - 50	5 - 50	1	2500.00
V - 200	50 - 200	1	2500.00
V - 1000	200 - 1000	1	2500.00
V - 100	10 - 100	1	2500.00
V - 500	100 - 1000	1	2500.00
V - 20	2 - 20	1	2500.00
V - 300	20 - 200	1	2500.00
V - 5000	50 - 5000	1	2500.00
V - 10000	500 - 1000	1	2500.00

Pipette Controller - Biopette

Features :

- Powerful pump with Double safety valve and Hydrophobic filter.
- Soft grip push buttons.
- Fill 25ml Pipette in just 4 seconds.
- Safety valve and hydrophobic filter provide double protection against fluid penetration.
- High and Low Speed, through varying Finger pressure.
- Better Control of speed. through varying Finger pressure.
- Low Battery Indicator.
- Cadmium free environment friendly batteries.
- Used Continuous for 8 Hours.
- Fully Battery Charging time 2-3 hrs.
- Pipette Volume Range from 0.1 - 200ml



Standard Accessories :

- Battery Charger
- Hydrophobic PTFE Membrane filter 0.45 μ l & 0.20 μ l each 1 No.
- Wall Mounting Stand with Screw

Cat No.	Description	Range	Price/Pc
PC-01	Pipette Controller	0.1 - 200ml	8500.00





LABORATORY PROTECTIVE GEARS



Eye Protection

Cat.No	Product Description	Unit Price
EP-01	<i>SAVIOUR Series 1 Safety goggle with Nylon frame punk design clear Polycarbonate lens. Ergonomically and sturdy design EN:166/EN:170</i>	70.00
EP-02	<i>SAVIOUR Series 2 Safety goggle with wrap around Polycarbonate lens. Ergonomically and sturdy design EN:166/EN:170</i>	70.00
EP-03	<i>3M 1621 Chemical Splash Goggles. Reduce exposure to a variety of eye hazards including impact, splash and dust. Features includes, indirect vents that helps protect against splash and a hard coated polycarbonate lens that reduces UVB and UV radiation from natural sunlight by 99%</i>	130.00
EP-04	<i>SAVIOUR EY-601 Safety Eyewear consists of a soft nylon frame and tough anti-scratch and anti-fog coated Polycarbonate lens for impact resistance along With high optical clarity and durability with adjustable arms EN:166</i>	150.00
EP-05	<i>3M VIRTUA IN-11-11850 Popular model unisex styled eyewear with a 3 piece polycarbonate construction having universal nose bridge and hard coated lens. ANSI Z 87.1 Approved</i>	90.00
EP-06	<i>3M 1709 IN Wide Frontal Coverage of eyes Lateral Protection. Clear lens for better clarity, Better fit due to reduced Temple Distance ANSI Z.87.1 Approved</i>	80.00



Hand Protection

Cat.No	Product Description	Unit Price
HP-01	LATEX Examination Gloves- Powdered Single use Pkt of 100	300.00
HP-02	LATEX Surgical Gloves Powder Free STERILE Box of 50Pairs	1100.00
HP-03	NITRILE Examination Gloves Powder Free Single use Pkt of 100	400.00
HP-04	NITRILE Gloves: Lakeland Chemical Resistance especially resistant to all hydrocarbons. Flock lined design for comfort Length 13" Palm size 9/1- EN:388/EN374	130.00
HP-05	LATEX/NEOPRENE Gloves: Ansell Protection against wide range of Chemicals. Double dipped for double Protection. Strong and durable with excellent	150.00
HP-06	NITRILE Gloves:Ansell Solvex Nitrile Glove is designed to deliver optimum result in wet and dry environments where chemical resistance is crucial Size 9/10/11 EN:388/EN374	135.00



Fire blankets:

FB-01

Fire Blankets:

Fire Blanket is made of fiber - glass fire retardant fabric in one meter by two meter length [1 x 2 mt]; edged sealed by fiber - glass tapes and corners provided with loops. The Fire Blanket is packed in a carrying bag with instructions printed.


Suggested Temperature Resistance : 200 - 300 °C





Respirator Protection

Cat.No	Product Description	Unit Price
RP-01	Face Mask 2Ply Non woven Disposable (Pkt of 100)	250.00
RP-02	Face Mask 3Ply Non woven Disposable (Pkt of 100)	300.00
RP-03	Respirator VENUS 420 stitched adjustable double elastic. No leakage. No Puncture on filter. Provides leak proof snug fit. Long Life, Skin friendly, does not deform in high temperature. Gives snug fit yet relieves strap pressure.	70.00
RP-04	Respirator VENUS 410SL Comfort folds Flat Style Respirator Ultrasonic welded to form composite filter without needle holes. One way exhalation vent valve.	40.00
RP-05	Saviour P2 Dust/Mist Respirator Disposable Respirators for P2-Class Protection.	20.00

Cat.No	Product Description	Unit Price
 <h2>Eye Wash Units</h2>		
EW-01	EyeWash bottle: Wall Mounted Stand is a solution for emergency flushing of eyes with water, eyesaling solution to wash away the hazardous chemicals. The bottle is compact and can be easily operated with one hand.	300.00
EW-02	Eye Wash Unit: Hand Operated Unit is made of Powder coated anti corrosive GI Pipes. The eye wash units are operated by a push valve. The Eyewash has a drain valve And is packed with all accessories separately	3700.00
EW-03	Eye Wash Unit: Foot Operated. Unit is made of Powder coated anti corrosive GI Pipes. The eye wash units are operated by a Foot Plate.The Eyewash has a drain valve And is packed with all accessories separately	4400.00

Terms of Sale

- ❖ Prices Mentioned in the Price List are the Maximum recommended Prices per Unit/ Piece in Indian Rupees.
 - ❖ Prices are Ex- Godown, Chennai, Packing, Forwarding Expenses extra at Cost. VAT or CST Extra as per Govt. rules.
 - ❖ Prices are Subject to change without Notice. Prices Prevailing at the time of actual delivery will be Valid.
 - ❖ Invoices are to be paid as per agreed Credit period without any deductions.
In Case of outstation Cheque / Drafts, The Bank Charges will be Customers account.
 - ❖ The Prices are for the Standard Pack for each item, it would be appreciated that Customer places orders for the items based on the Standard pack Size indicated and in multiples there of.
 - ❖ No goods should be returned to us without our Prior approval.
 - ❖ The Sales of the laboratory glassware products will be Subject to Chennai Jurisdiction for all the Parties and for all legal matters, the Indian Laws are applicable.
-
-



Product Range:

Autoclaves	BOD Incubators	Blood Storage Cabinets	BioSafety Cabinet
Carbon & Sulphur Determination Apparatus	Constant Temperature bath	Cryostat Bath	Bunsen burners
Environmental Chamber	Frozen/Semen Storage cabinets	Fume Exhaust Hoods	Deep Freezers (Chest & Upright)
Furnace (Grooved)	Furnace Tubular (1200°c)	Furnace High Temperature (1400°c)	Flocculator (Jar Test Apparatus)
Humidity and Temp. Cabinet	Hot Plates	Heating Mantles	Hot Air Oven
Ice Flaking Machine	Incinerator	Incubator	Humidifier
Incubator Shaker	Innoculatio hood	Kjeldahl (Digestion & Distillation)	Industrial Tray Drier
Low Temperature Bath	Mortuary Chamber	Muffle Furnace	Laminar Air Flow
Melting Point Apparatus	Orbital Shakers	Oil Baths	Magnetic Stirrers
Plant Growth Chamber	Rotary Shakers	Sieve Shakers	Oven(Industrial & Universal)
Sand baths	Tray Driers	Tissue Floatation bath	Soxhlet Extraction Units
VDRL Rotator	Vacuum Pump (Single./Double/Oil free)	Waterbaths	Vacuum Oven
Waterbath incubator shakers	Wrist action shakers	Water Still (wall & Table)	Waterbaths(Serological)



Manufactured & Marketed by:



Kavin Scientific Products

7A, Natarajan St., (Saidai Sambandan St.,)

Saidapet, Chennai-600015, INDIA

Tel: 044-43302975, 24714329

E-Mail: kavinscientific@yahoo.co.uk,

kavinscientific@gmail.com

Website: www.kavinscientific.in





THANK YOU